Performing the swab test for antibiotic residues (on premises)

**A self-instructional
guide**

**Agriculture
Canada**

©Minister of Supply and Services Canada 1980
Cat. No. A62-36/1980 ISBN: 0-662-10756-X
Printed 1980 200-6:80

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PERFORMING THE SWAB TEST (ON PREMISES) FOR ANTIBIOTIC RESIDUES

A self-instructional guide

The Canadian version of this instruction booklet was adapted with kind permission from the original developed by the United States Department of Agriculture, Food Safety and Quality Service, Meat and Poultry Inspection, Program Training Staff.

INTRODUCTION

A simple but highly accurate test for the detection of antibiotic residues in animal tissues has been developed. The test is called S.T.O.P. (Swab Test on Premises). The use of S.T.O.P. will reduce the holding time of many retained carcasses, reduce the time spent in packing and mailing specimens and will result in decreased mailing and laboratory costs.

Current instructions require that when the veterinarian in charge has reason to suspect that an animal's tissues may contain antibiotic residues, such as when injection lesions are found, the carcass is to be retained and appropriate samples sent to the laboratory for analysis. It often takes several days or more to receive the laboratory results. A high percentage of carcasses with visible injection lesions do not have detectable antibiotic residues in the tissues, meaning that after the time-consuming and expensive process of packing and mailing samples and a long holding period, the carcass is likely to be released.

By using S.T.O.P., the tissue may be tested in the slaughtering plant (on-premises) and the results obtained less than 24 hours after slaughter. If the S.T.O.P. results are negative, the carcass may be released immediately with a high degree of assurance that the tissues are free of antibiotic residues. Only S.T.O.P. positive cases will need to be submitted to the laboratory for confirmation and assay to determine the disposition.

S.T.O.P. is a biological test for the presence of inhibitors (antibiotics) in animal tissue. It is based on the principle that if the tissue contains an antibiotic residue, fluid from the tissue will inhibit the growth of a sensitive organism on a bacterial culture plate. Experience has shown that the kidney is the most likely organ to contain antibi-

otic residue, so kidney tissue is used for this test. The accuracy of S.T.O.P. has not been determined when used to test other tissues.

Cotton swabs saturated with kidney tissue fluid from the suspected carcass are placed on an agar gel culture plate. The agar gel plate has been preseeded with a special strain of *Bacillus subtilis* bacteria that is highly sensitive to all the common antibiotics. (NOTE: *Bacillus subtilis* is a harmless organism found in nature and is not dangerous to animal or human health.) The swabs and plate are incubated for 16-18 hours to allow growth of the organism. Then the plates are examined for zones of inhibited growth around the swabs.

Performing the S.T.O.P. is not difficult, does not require sterile technique, and takes only a few minutes of actual work time. The major steps in performing the S.T.O.P. are:

1. Preparing and setting up the agar gel plate;
2. Incubating the agar gel plate;
3. Reading, interpreting and recording test results; and
4. Taking and reporting action.

The performer must PREPARE tissue swabs and position them accurately in the agar gel plate.

The performer must INCUBATE the test plate for the required time and at the required temperature.

The performer must take proper ACTION based on the test results and REPORT them in a complete and timely manner.

This guide is a step-by-step procedure for performing the Antibiotic Swab Test on Premises (S.T.O.P.). If you follow the guide you will perform, interpret, and report the S.T.O.P. easily and accurately.

REMEMBER —

Use this guide each time you perform the S.T.O.P.

HOW TO USE THIS GUIDE

When you are assigned to test a retained carcass for the presence of antibiotic residues, take this guide to your work area.


Follow the instructions in this guide step-by-step until you have completed the test, recorded the results, taken required action, and submitted your report.

REMEMBER —


Use this guide each time you perform the antibiotic swab test on premises.

STEP 1

Collect approximately one-third of a kidney from the retained carcass and place the kidney tissue sample in a plastic bag.

Attach a held tag with the same number as the tag on the retained carcass to the kidney tissue sample bag.

SPECIAL NOTE: Make sure the carcass and its parts remain held until the S.T.O.P. is completed. Remember to obtain all carcass identification information [back tags, ear tags, owner identification, etc.].

Will you be on duty to read the S.T.O.P. results between 16 and 24 hours from now?

YES — Go to STEP 2.

NO — Retain the kidney tissue sample under refrigeration until you will be able to read the results 16-24 hours after starting the S.T.O.P. Then


CONTINUE

STEP 2

Locate a desk or table in the inspection office and clear a
60 × 60 cm work surface.

Obtain the following equipment from storage:

- One swab pack [cotton-tipped applicators]
 - Thumb forceps
 - S.T.O.P. measuring aid
 - S.T.O.P. report forms [with carbon]
 - Pencil
 - One sheet clean paper

Place the kidney tissue sample bag and equipment on the cleared work surface and examine the S.T.O.P. report form.

Are the results from a previous test recorded on the form?

YES — Go to STEP 3.

NO — Making an original and one copy, fill in the heading of the form. Then

CONTINUE.

STEP 3

Read the held tag attached to the kidney sample bag and write the held tag number in the first available block in the extreme left column of the S.T.O.P. report form.

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Select the code below that best describes the animal you are testing.

| | | |
|------------|-------------------|-------------------|
| Horse — 01 | Heifer — 05 | Goat — 09 |
| Bull — 02 | Calf — 06 | Market Swine — 10 |
| Steer — 03 | Mature Sheep — 07 | Boar or Stag — 11 |
| Cow — 04 | Lamb — 08 | Sow — 12 |

Write the species code you have selected in the block adjacent to the retained tag number.

Write today's date in the adjacent "In" block.

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
When you have entered the held tag number, species code, and date on the S.T.O.P. report form, draw a circle around the last two digits of the held tag number. Then

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
STEP 4


Pick up the swab pack and feel the swab shafts through the paper wrap to locate the approximate center of the shaft length. Then, grasp each side of the center and snap the plastic swab shafts.


Punch the broken end of the swab shafts through the paper wrap. Peel the paper wrap to expose the shafts. Pull one cotton-tipped swab from the pack and lay the swab pack aside.


Grasp the swab shaft just below the cotton wrapped tip and look through the plastic sample bag to locate one end or side of the kidney sample.


Jab the sharp broken end of the swab shaft through the plastic bag wall. Then use the swab shaft and the fingers of your other hand to enlarge the bag opening to about 12 mm in length.

Push the broken swab shaft about 12 to 19 mm into the kidney tissue. Then jab it back and forth several times to tear the kidney tissue and break down cells.


When the kidney tissue is torn and broken, withdraw the swab shaft, reverse it, and inset the cotton swab tip through the tear in the bag and into the torn kidney tissue.

Position the swab so that all of the cotton-wrapped portion is covered by kidney tissue and twirl the swab shaft to make sure there is good contact between the swab head and the kidney tissue. Leave the swab in place, and


CONTINUE.

STEP 5

Remove the other swab from the paper pack. Locate an area on the opposite end or other side of the kidney sample and

- puncture bag
- enlarge opening
- tear kidney tissue
- position and seat swab

When *both* swabs are seated in the kidney tissue, check the time and plan to leave the swabs in place at least 30 minutes but *not more than* 2 hours. This will ensure proper saturation of the swabs.

CONTINUE.


STEP 6

About 20 minutes before you are ready to remove the swabs from the kidney tissue, go to the refrigerated storage area and obtain the vial of N5 discs and one agar gel plate (with cover). *Be sure to reseal the plastic sleeve to prevent dehydration of other plates.*

Read the expiration date printed on the plate and compare that date to today's date.

Has the expiration date passed?

NO — Take the agar gel plate and the vial of N5 discs to the S.T.O.P. work area and
go to STEP 7.

YES — Discard the outdated agar gel plate into the trash. Examine other plates in storage and discard all outdated plates.

Are any plates that are NOT outdated available?

YES — Take a current plate and the vial of N5 discs to the S.T.O.P. work area and
go to STEP 7.

NO — When all outdated plates are discarded, return to the S.T.O.P. work area and
go to STEP 22.

STEP 7

SPECIAL NOTE: Each time you set up the S.T.O.P., you must first determine that the agar gel plate has been properly refrigerated during storage. The next two steps will help you do this.


Wash, rinse, and dry your hands.

Remove the cover from the agar gel plate and hold the plate at an angle so you can observe the surface of the agar gel layer.

Is the surface of the agar gel dry and shiny [about like very firm gelatin]?

YES — Go to STEP 8.


NO — CONTINUE.

Is there free water on the surface of the agar gel or does the surface appear wrinkled?

YES — The plate has been frozen during storage and is unsuitable. Discard the plate and **go back to STEP 6.**


NO — CONTINUE.

STEP 8

Turn the agar gel plate over and hold it with the bottom of the plate toward you so you can look through the plate at the agar gel layer.


SPECIAL NOTE: Lighting conditions vary so you will have to experiment to hold the plate in a way that gives you the best possible view of the agar gel layer. One good method is to hold the plate up to a shaded desk lamp so the light shines through the plate but not directly into your eyes. If the plate is still cold you may have to wipe condensation from the bottom of the plate to see the gel layer clearly.

Is the agar gel layer a uniform, slightly opaque, light amber color, with no visible bacterial colony growth?

YES — Consider the plate acceptable and
go to STEP 9.


NO — CONTINUE.

Are there many white pinpoint colonies uniformly distributed across the agar gel layer?

YES — The plate has become too warm during storage and cannot be used. Discard the plate and
go back to STEP 6.

NO — CONTINUE.

There may be one or more fairly large colonies [possibly pigmented] growing on the surface of the agar gel layer. These are the growth of some stray organism that grows at low temperature [psychrophile].

Are there several of these large colonies or is there one colony near the centre of the plate?

YES — These colonies may interfere with the test. Discard the plate and
go back to STEP 6.


NO — One or two colonies near the edge of the plate will not interfere with the test. Plan to avoid these colonies when placing the N5 disc and swabs on the plate.
CONTINUE.

STEP 9


Locate and read the last two digits of the kidney tissue held tag number.

Write this two-digit number in the space provided on the agar gel plate label.

Check to make sure the plate number matches the circled number on the S.T.O.P. report form.


Replace the cover on the agar gel plate and

CONTINUE.

STEP 10


Position the clean sheet of paper on the work surface and place the S.T.O.P. measuring aid on the paper.

With the cover on top, position the agar gel plate on the measuring aid so the plate is within the circle.


Open the glass vial and remove the N5 disc dispenser


SPECIAL NOTE: Leave the small dehydrator capsule in the glass vial. This prevents the discs from picking up moisture during storage.


Being careful not to touch the N5 discs with your fingers, remove the plastic cap from the N5 disc dispenser and set the cap aside.


Lift the cover from the agar gel plate and place it open side up beside the plate.

Hold the N5 disc dispenser over the open cover and dispense one N5 disc into the cover by pressing the plastic plunger.


Replace the plastic cover on the dispenser, put it back into the glass vial, and set the vial aside.


Using the thumb forceps, pick up the N5 disc by its edges, and hold the disc just above the surface of the agar gel layer directly over the center measuring aid dot. Release finger pressure on the forceps to drop the N5 disc flat onto the approximate center of the agar gel layer.

SPECIAL NOTE: *Leave the N5 disc where it drops.* DO NOT attempt to reposition the disc nearer to the center of the agar gel plate. DO NOT touch the N5 disc with your fingers.


Being careful *not* to press hard enough to break the surface of the agar gel layer, *lightly* touch the disc with the forceps to assure uniform contact.


Is the N5 disc seated in the approximate center of the agar gel plate [at least partly covering the center dot]?

YES — Replace the cover on the agar gel plate and
go to STEP 11.


NO — Discard the agar gel plate and
go back to STEP 6.


STEP 11

Check the time. When the swabs have been in the kidney tissue at least 30 minutes,


Wash, rinse and dry your hands. Then, pull one swab from the kidney tissue, hold it with the cotton tip down and grasp the plastic shaft with your *left* thumb and forefinger just above the cotton. DO NOT touch the cotton tip itself.


Apply bending pressure with the fingers of the other hand to break the swab shaft short.

Transfer the short swab to your right hand by grasping the broken tip of the shaft with your thumb and forefinger and lift the cover with the other hand.


Hold the short swab just above the agar gel layer, align it with the top measuring aid line and gently drop it onto the agar gel layer.


Being careful not to break the agar gel surface, *lightly* press the short swab with your finger tip where the swab head joins the shaft to assure uniform contact of the swab head with the agar gel layer.

SPECIAL NOTE: The thumb forceps are NOT used to seat the swab because the forceps may be contaminated with neomycin from handling the N5 disc.


When the first swab has been seated on the agar gel plate,


rotate the agar gel plate one-half turn so the swab shaft is aligned with the lower measuring aid line. Replace the cover. Then

CONTINUE.


STEP 12


Pull the second swab from the kidney tissue, break the shaft short, and lift the cover.

Again, using the upper measuring aid line as a guide, position the second short swab halfway between the N5 disc and the edge of the plate on the opposite side from the first swab.


Gently seat the second swab with your fingertip.

When both swabs are positioned and seated, *replace the cover* and

CONTINUE.

STEP 13

Check to make sure the agar gel plate is labeled with the last two digits of the held tag number. Then, place the test plate in a stabilized incubator [27-29°C] with the cover on top.

Secure the incubator as required.

Note the time now and on the S.T.O.P. report form, write the time in the "Time In" block next to the date.

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Plan to incubate the test plate overnight.

SPECIAL NOTE: The ideal incubation time for S.T.O.P. is about 16-18 hours. The results may be difficult to read before 16 hours and prolonged incubation [more than 24 hours] may result in a false negative test.

Return the kidney tissue sample and the vial of N5 discs to *secure refrigerated storage*.

Return other equipment to supply storage.

When incubation is complete,

CONTINUE.


STEP 14

Select the following from supply storage:

- the partially completed S.T.O.P. report form
- S.T.O.P. measuring aid
- pencil

Remove the test plate from the incubator and place it on the table with the *cover down* and the agar gel plate up.


SPECIAL NOTE: If you are unable to read the results immediately after the 16-18 hour incubation, or need assistance in interpreting the results, you may place the test plate in a refrigerator to stop further growth and read the results later.

Locate the "Incubation" column on the S.T.O.P. report form and on the line where the circled held tag digits match the plate number, enter today's date and time the test plate was taken out of the incubator in the "Out" blocks.

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
Tap the bottom of the agar gel plate with your fingertip until the swabs are dislodged from the agar gel surface and drop into the inverted cover.

Lift the agar gel plate out of the inverted cover. Hold the plate with the bottom of the plate toward you so you can look through the plate at the agar gel layer.

SPECIAL NOTE: Position the plate so you can get the best possible view of the agar gel layer with your lighting conditions.


Examine the agar gel layer for the presence of bacterial colony growth.

Are there many tiny white colonies uniformly distributed in at least part of the agar gel layer?

YES — This is growth of the seeded organism.
Go to STEP 15.

NO — CONTINUE.

Except for discolored areas where the swabs were in contact with the agar gel, is the agar gel layer a uniform, slightly opaque, light amber color [as it was before incubation] with no bacterial colonies present?

YES — Go to STEP 24.


NO — CONTINUE.

There may be a few fairly large randomly distributed colonies [possibly pigmented] on the surface of the agar gel layer. These are the growth of some stray organism.


In the absence of many tiny colonies uniformly distributed in the agar gel layer, consider these larger surface colonies the same as no culture growth and

go to STEP 24.

STEP 15


Disregard any large or pigmented surface colonies regardless of their location. These are from the growth of some stray organism and have no significance. Consider only the *tiny closely spaced uniformly distributed white colonies* as evidence of growth of the seeded organism.


Examine the area immediately around the N5 disc.

Is there a clear zone [inhibition zone] around the N5 disc that is free of colonies of the seeded organism?

YES — Go to STEP 16.

NO — CONTINUE.


Locate the "N5 Disc Zone" column on the S.T.O.P. report form and write "None" in the appropriate block. Then

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STEP 16

Pick up the S.T.O.P. measuring aid, position it against the bottom of the agar gel plate with the millimeter scale across the N5 disc parallel to the swab impressions, and measure the diameter of the clear zone.

Locate the "N5 Disc Zone" column on the S.T.O.P. report form and record the measured diameter in the appropriate block.

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
STEP 17

Examine one of the swab impressions in the agar gel layer to see if there is a clear zone free of colonies of the seeded organism surrounding the area where the swab head [cotton tip] was in contact with the agar gel.


SPECIAL NOTE: The swab head impression area may be slightly colored due to diffusion of kidney tissue pigments from the swab head into the agar gel layer, so look carefully for growth. Also, pressure from the swab may prevent colonies from developing normally. DO NOT consider lack of growth where the swab actually contacted the agar gel as a clear zone.

Is there a clear zone surrounding the first swab impression?

YES — Go to STEP 18.

NO — CONTINUE.


Locate the appropriate double block in the "Swab Zones" column on the S.T.O.P. report form and write "None" in the upper half of this double block.

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
STEP 18

Using the S.T.O.P. measuring aid, determine the length and width of the clear zone in millimeters.

SPECIAL NOTE: If the clear zone around the swab impression overlaps the clear zone around the N5 disc, measure the width from the center of the swab impression to the edge of the zone and double this measurement.

Locate the appropriate double block in the "Swab Zones" column on the S.T.O.P. report form and write the length and width measurements in the *upper half* of this double block. Record these as "length measurement X width measurement" in millimeters.

When the first swab zone dimensions are recorded,

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STEP 19

Examine the other swab impression area.

Is there a clear zone surrounding the second swab impression?

YES — Measure the length and width of the clear zone in millimeters and record the measurements in the *lower half* of the appropriate "Swab Zones" block.
Then

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NO — Write "None" in the *lower half* of the appropriate "Swab Zones" block. Then

CONTINUE.

STEP 20

Examine the double block entries in the "Swab Zones" column on the S.T.O.P. report form.

Is there a length × width swab head clear zone measurement in either [or both] of the blocks?

| N5 Disc Zone (MM) | Surface du Disque NB (MM) | Swab Zones L x W MM Surfaces écouvillonées Long x Large (MM) | Test Results Résultats |
|-------------------|---------------------------|---|------------------------|
| 18 | | <u>NONE</u> | |
| 20 | | <u>NONE</u> | <u>N</u> |
| 19 | | <u>9x5</u> <u>11x7</u> | <u>Pt+</u> |
| | | <u>NONE</u> | |
| | | <u>8x6</u> | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

| N5 DISC Zone (MM) Surface du Disque N5 (MM) | Swab Zones L x W MM Surfaces écouvelloées Long x Large (MM) | Test Résultats |
|---|---|-------------------|
| 18 | <u>NONE</u> | N |
| 20 | <u>NONE</u> <u>9x5</u> <u>11x7</u> | P+ |
| 19 | <u>6x4</u> <u>8x6</u> | |
| | | |
| | | |
| | | |
| | | |
| | | |

YES — Consider the test result to be PRESUMPTIVE POSITIVE. Write "P+" in the "Test Results" block for this tissue sample. Then

go to STEP 25.

NO — Go to STEP 21.

STEP 21

You have determined that the S.T.O.P. test result was NEGATIVE [there was no growth inhibition zone in either swab impression area].

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Calculate the number of hours the test plate was incubated.
Was the test plate incubated *more than* 24 hours?

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NO — Write "N" in the "Test Results" block for this tissue sample.

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Recall the reason the S.T.O.P. was performed on this carcass [injection lesion, mastitis, pneumonia, etc.].

Write this reason in the appropriate block in the "Remarks" column on the S.T.O.P. report form.

Unless the carcass was retained for reasons other than the S.T.O.P., release the carcass and its parts.

SPECIAL NOTE: Be sure to require trimming of any injection lesions that may be present.

Return the S.T.O.P. measuring aid to storage and discard the test plate. Then

go to STEP 26.

STEP 22

You have found that no current-dated agar gel plates are available. You will be unable to perform the S.T.O.P. until a new supply of plates arrives.

Contact your regional office to order a new supply of agar gel plates.

On the S.T.O.P. report form, locate the "N5 Disc Zone" and "Swab Zones" columns and write "No Current Plates" across both columns on the line for this tissue sample to indicate that you are unable to perform the S.T.O.P.

When you have recorded that you are unable to test this tissue sample, remove and discard the two swabs and

go to STEP 25.

STEP 23

Prolonged incubation [more than 24 hours] reduces the accuracy of the test. Inhibition activity in the swab may dissipate, allowing growth of the organism, thus giving a false negative test result.

Locate the "Test Results" column on the S.T.O.P. report form and write "I" in the "Test Results" block for this sample to indicate the inconclusive S.T.O.P. results.

| In | N5 DISC Zone (MM) Surface du Disque N5 (MM) | Swab Zones L x W MM Surfaces écouvillonnées Long x Large (MM) | Test Results Résultats | Tissue Samples sent To: Echantillons envoyés à |
|---------------|---|--|---------------------------|---|
| Time Heure | | | | |
| 13:15 | 18 | NONE | | |
| 9:00 | | NONE | N | - |
| 13:30 | | 9x5 | | |
| 9:00 | 20 | 11x7 | P+ | ADRI(E) |
| 15:30 | | NONE | | |
| 11:00 | 19 | NONE | I | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Decide if you wish to rerun the S.T.O.P. or submit the carcass tissues to the laboratory.

IF you have decided to rerun the S.T.O.P.,

go back to STEP 1.

IF you have decided to submit the carcass tissues to the laboratory,

go to STEP 25.

STEP 24

You have determined that there was no growth of the seeded organism.

Locate the columns "N5 Disc Zone" and "Swab Zones" on the S.T.O.P. report form. Write "No culture growth" across both columns on the line where the circled retained tag digits match the plate number.

Culture growth failure is a rare occurrence. Contact:

Chief, Residues
Meat Hygiene Directorate
Health of Animals Branch
Sir William Logan Building
580 Booth Street
Ottawa, Ontario

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STEP 25

You have determined that

- no current-dated agar gel plates were available
 - or
 - the test plate was incubated more than 24 hours
 - or
 - there was no culture growth
 - or
 - the test result was PRESUMPTIVE POSITIVE [P+].

Maintain retention of the carcass and follow instruction in Directive of "date" to collect, prepare, and pack tissues from the retained carcass for shipment to the laboratory.

Complete Form HA89 to accompany the tissues to the laboratory. In addition to the regular entries on the HA89, make one of the following entries in the Inspector's Remarks section.

If the test plate was incubated more than 24 hours, write "S.T.O.P. incubated more than 24 hours."

If there were no current-dated agar gel plates available, write "No Current S.T.O.P. plates available."

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When completing HA89, insert held tag number, reason for test and N if negative or P+ if positive.

Ship the packed samples to the laboratory servicing your area.

Locate the "Tissue Samples Sent To" column on the S.T.O.P. report form. Write the name of the laboratory in the appropriate block in this column.

When the laboratory city has been entered,

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STEP 26

Is the S.T.O.P. report form full with no additional lines available for use?

YES — Go to STEP 27.

NO — Continue to use the form until it is full or until the first day of the next fiscal quarter [Jan. 1, Apr. 1, Jul. 1, Oct. 1], whichever occurs first.

CONTINUE.

STEP 27

Mail the original S.T.O.P. report form to:

Chief, Residues
Meat Hygiene Directorate
Health of Animals Branch
Sir William Logan Building
580 Booth Street
Ottawa, Ontario

File the duplicate copy in your file.

SPECIAL NOTE: Check your inventory of S.T.O.P. supplies and order any needed supplies through your regional office.

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